## **GF-P Series** Pharmacy Balances

The GF-P Series balances have met the rigorous standards set forth by the National Conference on Weights and Measures Handbook 44, for legal-for-trade use. The GF-P balances are ideal for pharmacy and retail applications or any situation that requires legal-for-trade precision balances.

- Ideal for Pill Counting & Compounding
- Counting feature for prescription filling
- Meets Handbook 44 Requirements for a "Class A" Prescription Balance
- Legal-for-Trade Prescription Balance
- IP-54 splash proof operator keyboard protects against dust and spills, stainless steel weigh pan
- GLP/GMP/ISO compliance
- One-Second Response Time





## **GF-P SERIES SPECIFICATIONS**

Model		GF-200P	GF-300P
Weighing capacity		210 g	310 g
Readability		0.001 g (d) value 0.01 g (e) value (NTEP Class II)	
Repeatability		0.001g (Standard Deviation)	
Other units of measure		Decimal ounce, Troy ounce, Pennyweight, Grain unit, Pound, Pound/Ounce	
Linearity		±0.002 g	
Stabilization time (typical at FAST)		Approx. 1 second	
Sensitivity drift (10°C to 30°C/50°F-80°F)		±2ppm/°C (When automatic self calibration is not used)	
Operating temperature		5°C to 40°C (41°F to 104°F), 85%RH or less (No condensation)	
Sensing method		Super Hybrid Sensor (SHS)	
Display type		Vacuum Fluorescent Display (VFD)	
Display refresh rate		5 times/second or 10 times/second	
Counting mode	Minimum unit mass	0.001 g	
	Number of samples	10, 25, 50 or 100 pieces	
Percent mode	Minimum 100% reference mass	0.100 g	
	Minimum 100% display	0.01%, 0.1%, 1% (Depends on reference mass stored)	
Standard serial I/F		RS-232C Interface with Windows Communication Tools (WinCT: included as CD-ROM)	
Weighing pan		128 x 128 mm (5.04 x 5.04 inch)	
Calibration		External Calibration	
External dimensions		210(W) x 317(D) x 86(H) mm (8.27 x 12.48 x 3.39 inch)	
AC adapter		Confirm that the adapter type is correct for the local voltage and power receptacle	
Power consumption		Approx. 12W (supplied by AC adapter)	
Net weight		Approx. 8.5 lbs	

## **Options**

- GX-02\* Uni-directional Quick USB interface with cable
- GX-04\* Comparator Output with a Buzzer/RS-232C/ **Current Loop Output**
- GX-06\* Analog Output/Current Loop Output
- GX-10 Tall Breeze Break for GF-200P & GF-300P
- GX-12 Animal weighing pan for GF-300P
- GX-13 Density Determination Kit for GF-200P & GF-300P \* If selected, only one of option GX-02, GX-04 or GX-06 can be installed (factory installed option).

## Accessories

- AD-1192
- KO:1131
- AD-8121B
- AD-8920A
- **Rechargeable Battery Unit** • AD-1682
  - (Maximum operating time: 8 hours) RS-232C Cable (25P-9P)

**Dot Matrix Compact Printer** 

- KO:WW9-25
- SW-128
- Foot Switch • AX-073003691-S Dust Cover for GF-200 & GF-300P

**Compact Printer** Cable for AD-1192 (3 ft.)

**Remote Display** 



1756 Automation Parkway San Jose, CA 95131 (800) 726-3364 or (408) 263-5333 Email: scales@andweighing.com www.andweighing.com



Specifications subject to change without notice.